

IL17A antibody [MT504] (Biotin)

Cat. No. GTX02966-02

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Application	FACS, ELISA, Sandwich ELISA
Reactivity	Human, Sheep, Bovine, Dog, Pig, Monkey, Horse

Package
250 µg

APPLICATION

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
FACS	Assay dependent
ELISA	Assay dependent
Sandwich ELISA	Assay dependent
Note : Capture: GTX02964 / GTX02965, Detection: GTX02966-02.	

Not tested in other applications.

Calculated MW	18 kDa. (Note)
Product Note	This antibody is able to detect native and recombinant human IL-17A.

PROPERTIES

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant human IL-17A
Purification	Protein G purified From tissue culture supernatant
Conjugation	Biotin Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTx02966-02 ELISA Image

Sandwich ELISA analysis of monkey IL-17A protein using GTX02964 IL17A antibody [MT241] as coating antibody and GTX02966-02 IL17A antibody [MT504] (Biotin) as detecting antibody.

Substrate : pNPP


GTx02966-02 ELISA Image

Sandwich ELISA analysis of human IL-17A protein using GTX02965 IL17A antibody [MT44.6] as coating antibody and GTX02966-02 IL17A antibody [MT504] (Biotin) as detecting antibody.

Substrate : pNPP



For full product information, images and publications, please visit our [website](#).