

IL21 antibody [MT216G]

Cat. No. GTX02978

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Application	ELISA, Sandwich ELISA
Reactivity	Human, Monkey

Package
250 µg

APPLICATION

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	Assay dependent
Sandwich ELISA	Assay dependent

Note : Capture: GTX02978, Detection: GTX02977-02.

Not tested in other applications.

Calculated MW	18 kDa. (Note)
Product Note	This antibody is able to detect native and recombinant human IL-21.

PROPERTIES

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant human IL-21
Purification	Protein G purified From tissue culture supernatant
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).


Date 2024 / 05 / 20 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTx02978 ELISA Image

Sandwich ELISA analysis of human IL-21 protein using GTx02978 IL21 antibody [MT216G] as coating antibody and GTx02977-02 IL21 antibody [MT21.3m] (Biotin) as detecting antibody.

Substrate : pNPP



For full product information, images and publications, please visit our [website](#).