

IL22 antibody [MT12A3]

Cat. No. GTX02979

Host	Rat
Clonality	Monoclonal
Isotype	IgG2a
Applications	ELISA, Sandwich ELISA
Reactivity	Human

Package
250 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	Assay dependent
Sandwich ELISA	Assay dependent

Note : Capture: GTX02979, Detection: GTX02982-02.

Not tested in other applications.

Product Note This antibody is able to detect native and recombinant human IL-22.

Properties


Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant human IL-22
Purification	Protein G purified From tissue culture supernatant
Conjugation	Unconjugated

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX02979 ELISA Image

Sandwich ELISA analysis of human IL-22 protein using GTX02979 IL22 antibody [MT12A3] as coating antibody and GTX02982-02 IL22 antibody [MT7B27] (Biotin) as detecting antibody.



For full product information, images and publications, please visit our [website](#).