

IL27 antibody [MT27 + MT361]

Cat. No. GTX02984

Host	Mouse
Clonality	Multiclonal
Isotype	IgG1, IgG2b
Applications	ELISA, Sandwich ELISA
Reactivity	Human

Package
250 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	Assay dependent
Sandwich ELISA	Assay dependent

Note : Capture: GTX02984 (MT27 + MT361 for IL27) / GTX02920 (MT23A14 for EBI3), Detection: GTX02919-02.

Not tested in other applications.

Product Note

This antibody is able to detect native and recombinant human IL-27. Clone [MT27] binds to IL-27, and clone [MT361] binds to p28, which is part of the IL-27 heterodimer.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant human IL-27
Purification	Protein G purified From tissue culture supernatant
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTX02984 ELISA Image

Sandwich ELISA analysis of human IL-27 protein using GTX02984 IL27 antibody [MT27 + MT361] as coating antibody and GTX02919-02 EBI3 antibody [MT140] (Biotin) as detecting antibody.



For full product information, images and publications, please visit our [website](#).