

## IL4 antibody [bIL4-I]

Cat. No. GTX02990

<b>Host</b>	Mouse
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	IgG2a
<b>Applications</b>	ELISA, Sandwich ELISA
<b>Reactivity</b>	Sheep, Bovine

Package  
250 µg

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	Assay dependent
Sandwich ELISA	Assay dependent

**Note : Capture: GTX02990, Detection: GTX02991-02.**

Not tested in other applications.

**Product Note** This antibody is able to detect native and recombinant bovine IL-4.

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	PBS
<b>Preservative</b>	0.02% Sodium azide
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
<b>Concentration</b>	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	Recombinant bovine IL-4
<b>Purification</b>	Protein G purified From tissue culture supernatant
<b>Conjugation</b>	Unconjugated


**Note** For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

## DATA IMAGES

**GTX02990 ELISA Image**

Sandwich ELISA analysis of bovine IL-4 protein using GTX02990 IL4 antibody [bIL4-I] as coating antibody and GTX02991-02 IL4 antibody [bIL4-II] (Biotin) as detecting antibody.

Substrate : TMB



For full product information, images and publications, please visit our [website](#).