

## IL5 antibody [5A10] (Biotin)

**Cat. No. GTX02995-02**

<b>Host</b>	Rat
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	IgG2a
<b>Applications</b>	ELISA
<b>Reactivity</b>	Human, Monkey

**Package**  
250 µg

## Applications

**Application Note**

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	Assay dependent

Not tested in other applications.

**Product Note** This antibody is able to detect native and recombinant human IL-5.

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	PBS
<b>Preservative</b>	0.02% Sodium azide
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
<b>Concentration</b>	1 mg/ml (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	Recombinant human IL-5
<b>Purification</b>	Protein G purified From tissue culture supernatant
<b>Conjugation</b>	Biotin Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.


**Note** For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

## DATA IMAGES

**GTX02995-02 ELISA Image**

Sandwich ELISA analysis of human IL-5 protein using GTX02997 IL5 antibody [TRFK5] as coating antibody and GTX02995-02 IL5 antibody [5A10] (Biotin) as detecting antibody.



For full product information, images and publications, please visit our [website](#).