

IL12 p70 antibody [MT704] (Biotin)

Cat. No. GTX03013-02

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Applications	ELISA, Sandwich ELISA
Reactivity	Human, Monkey

Package
250 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	Assay dependent
Sandwich ELISA	Assay dependent

Note : Capture: GTX02957, Detection: GTX03013-02 (MT704 for p70) / GTX02956-02 (MT618 for p40) / GTX02983-02.

Not tested in other applications.

Product Note

This antibody is specific for IL-12 p70 (p35 + p40). In combination with MT86/221 (p40 specific), it detects IL-12 p70 heterodimer (p35 + p40) by sandwich ELISA.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant human IL-12 p70
Purification	Protein G purified From tissue culture supernatant
Conjugation	Biotin Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

**GTX03013-02 ELISA Image**

Sandwich ELISA analysis of human IL-12 (p70) protein using GTX02957 IL12B / IL12 p40 antibody [MT86 + MT221] as coating antibody and GTX03013-02 IL12 p70 antibody [MT704] (Biotin) as detecting antibody.



For full product information, images and publications, please visit our [website](#).