

Interferon alpha antibody [MT1]

Cat. No. GTX03014

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	ELISA, Sandwich ELISA
Reactivity	Human

Package
250 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	Assay dependent
Sandwich ELISA	Assay dependent

Note : Capture: GTX03014, Detection: GTX03016-02.

Not tested in other applications.

Product Note This antibody is able to detect native and recombinant human IFN- α subtypes 2a, 2b, and 2c.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Human recombinant Interferon- α 2c
Purification	Protein G purified From tissue culture supernatant
Conjugation	Unconjugated


Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX03014 ELISA Image

Sandwich ELISA analysis of human IFN- α type 2 protein using GTX03014 Interferon alpha antibody [MT1] as coating antibody and GTX03016-02 Interferon alpha antibody [MT2] (Biotin) as detecting antibody.



For full product information, images and publications, please visit our [website](#).