

Mouse Anti-Human IgM antibody [MT11 + MT12]

Cat. No. GTX03024

| | |
|---------------------|---------------------------|
| Host | Mouse |
| Clonality | Multiclonal |
| Isotype | IgG1, IgG2b |
| Applications | WB, ELISA, Sandwich ELISA |
| Reactivity | Human, Monkey |

Package
250 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB | Assay dependent |
| ELISA | Assay dependent |
| Sandwich ELISA | Assay dependent |

Note : Capture: GTX03024, Detection: GTX03025-02.

Not tested in other applications.

Properties

| | |
|----------------------|--|
| Form | Liquid |
| Buffer | PBS |
| Preservative | 0.02% Sodium azide |
| Storage | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| Concentration | 0.5 mg/ml (Please refer to the vial label for the specific concentration.) |
| Immunogen | Human IgM |
| Purification | Protein G purified From tissue culture supernatant |
| Conjugation | Unconjugated |


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX03024 ELISA Image

Sandwich ELISA analysis of human IgM protein using GTX03024 Human IgM antibody [MT11 + MT12] as coating antibody and GTX03025-02 Human IgM antibody [MT22] (Biotin) as detecting antibody.



For full product information, images and publications, please visit our [website](https://www.genetex.com).