

Rat anti-Mouse IgE antibody [MT18E]

Cat. No. GTX03028

Host	Rat
Clonality	Monoclonal
Isotype	IgG2a
Applications	WB, ELISA, Sandwich ELISA
Reactivity	Mouse

Package
250 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ELISA	Assay dependent
Sandwich ELISA	Assay dependent

Note : Capture: GTX03028 / GTX03030, Detection: GTX03029-02.

Not tested in other applications.

Product Note The antibody is specific for mouse IgEa and does not detect allotype b, mouse IgEb.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Mouse IgE
Purification	Protein G purified From tissue culture supernatant
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTX03028 ELISA Image

Sandwich ELISA analysis of mouse IgEa (IgE allotype a) protein using GTX03028 Rat anti-Mouse IgE antibody [MT18E] as coating antibody and GTX03029-02 Rat anti-Mouse IgE antibody [MT44E] (Biotin) as detecting antibody.

Substrate : TMB



For full product information, images and publications, please visit our [website](#).