

Rat anti-Mouse IgM antibody [MT6A3]

Cat. No. GTX03031

Host	Rat
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ELISA, Sandwich ELISA
Reactivity	Mouse

Package
250 µg

Applications
Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ELISA	Assay dependent
Sandwich ELISA	Assay dependent

Note : Capture: GTX03031, Detection: GTX03032-02.

Not tested in other applications.

Product Note This antibody is able to detect native mouse IgM.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Mouse IgM
Purification	Protein G purified From tissue culture supernatant
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

**GTX03031 ELISA Image**

Sandwich ELISA analysis of mouse IgM protein using GTX03031 Rat anti-Mouse IgM antibody [MT6A3] as coating antibody and GTX03032-02 Rat anti-Mouse IgM antibody [MT9A2] (Biotin) as detecting antibody.



For full product information, images and publications, please visit our [website](#).