

## CXCL8 / IL8 antibody [MT8H6]

Cat. No. GTX03037

<b>Host</b>	Mouse
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	IgG1
<b>Applications</b>	ELISA, Sandwich ELISA
<b>Reactivity</b>	Human, Bovine, Dog, Cynomolgus monkey

Package  
250 µg

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	Assay dependent
Sandwich ELISA	Assay dependent

**Note : Capture: GTX03037, Detection: GTX02917-02 / GTX02918-02.**

Not tested in other applications.

**Product Note** This antibody is able to detect native and recombinant bovine IL-8.

## Properties


<b>Form</b>	Liquid
<b>Buffer</b>	PBS
<b>Preservative</b>	0.02% Sodium azide
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
<b>Concentration</b>	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	Recombinant human IL-8
<b>Purification</b>	Protein G purified From tissue culture supernatant
<b>Conjugation</b>	Unconjugated

**Note** For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

**DATA IMAGES**

**GTX03037 ELISA Image**

Sandwich ELISA analysis of human IL-8 (CXCL8) protein using GTX03037 CXCL8 / IL8 antibody [MT8H6] as coating antibody and GTX02918-02 CXCL8 / IL8 antibody [MT8F19] (Biotin) as detecting antibody.


**GTX03037 ELISA Image**

Sandwich ELISA analysis of bovine IL-8 (CXCL8) protein using GTX03037 CXCL8 / IL8 antibody [MT8H6] as coating antibody and GTX02917-02 CXCL8 / IL8 antibody [26E5] (Biotin) as detecting antibody.



For full product information, images and publications, please visit our [website](#).