

Bovine IgG ELISA pair [MT134/MT391]

Cat. No. GTX03040
Applications ELISA, Sandwich ELISA

Reactivity Bovine

Package

1 pair

PRODUCT

Product Content

Cat No.	Type	Host	Conjugation	Package
GTX03033	Capture antibody	Rat	Unconjugated	250 µg
GTX03034-02	Detection antibody	Rat	Biotin	250 µg

Properties

Form Liquid

Storage


Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTX03040 ELISA Image

Sandwich ELISA analysis of bovine IgG protein using GTX03033 Rat anti-Bovine IgG antibody [MT134] as coating antibody and GTX03034-02 Rat anti-Bovine IgG antibody [MT391] (Biotin) as detecting antibody. Substrate : TMB



For full product information, images and publications, please visit our [website](#).