

Mouse IgM ELISA pair [MT6A3/MT9A2]

Cat. No. GTX03044

Applications ELISA, Sandwich ELISA

Reactivity Mouse

Package

1 pair

PRODUCT

Product Content

Cat No.	Type	Host	Conjugation	Package
GTX03031	Capture antibody	Rat	Unconjugated	250 µg
GTX03032-02	Detection antibody	Rat	Biotin	250 µg

Properties


Form Liquid

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTX03044 ELISA Image

Sandwich ELISA analysis of mouse IgM protein using GTX03031 Rat anti-Mouse IgM antibody [MT6A3] as coating antibody and GTX03032-02 Rat anti-Mouse IgM antibody [MT9A2] (Biotin) as detecting antibody.



For full product information, images and publications, please visit our [website](#).