

## Apolipoprotein E ELISA pair [E276/E887]

Cat. No. GTX03050

Application ELISA, Sandwich ELISA

Reactivity Human

Package

1 pair

## PRODUCT

## Product Content

Cat No.	Type	Host	Conjugation	Package
GTX02904	Capture antibody	Mouse	Unconjugated	250 µg
GTX02905-02	Detection antibody	Mouse	Biotin	250 µg

## PROPERTIES

Form Liquid

## Storage


Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

## DATA IMAGES



## GTX03050 ELISA Image

Sandwich ELISA analysis of human apoE protein using GTX02904 Apolipoprotein E antibody [E276] as coating antibody and GTX02905-02 Apolipoprotein E antibody [E887] (Biotin) as detecting antibody.



For full product information, images and publications, please visit our [website](#).