

# Human IgA (Fc) ELISA pair [MT57/MT20]

**Cat. No. GTX03062**

**Application** ELISA, Sandwich ELISA

**Reactivity** Human

**Package**

1 pair

## PRODUCT

### Product Content

Cat No.	Type	Host	Conjugation	Package
GTX03018-02	Detection antibody	Rat	Biotin	250 µg
GTX03019	Capture antibody	Rat	Unconjugated	250 µg

## PROPERTIES

**Form** Liquid

### Storage


Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

### Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

## DATA IMAGES



### GTX03062 ELISA Image

Sandwich ELISA analysis of human IgA protein using GTX03019 Human IgA antibody [MT57] as coating antibody and GTX03018-02 Human IgA antibody [MT20] (Biotin) as detecting antibody.

Substrate : pNPP



For full product information, images and publications, please visit our [website](#).