

# IL1 alpha ELISA pair [MT641/642/MT513]

**Cat. No. GTX03066**

**Applications** ELISA, Sandwich ELISA

**Reactivity** Human

**Package**

1 pair

## PRODUCT

### Product Content

Cat No.	Type	Host	Conjugation	Package
GTX02967-02	Detection antibody	Mouse	Biotin	250 µg
GTX02968	Capture antibody	Mouse	Unconjugated	250 µg

## Properties


**Form** Liquid

**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

**Note** For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

## DATA IMAGES



### GTX03066 ELISA Image

Sandwich ELISA analysis of human IL-1α protein using GTX02968 IL1 alpha antibody [MT641 + MT642] as coating antibody and GTX02967-02 IL1 alpha antibody [MT513] (Biotin) as detecting antibody.



For full product information, images and publications, please visit our [website](#).