

IL-12 p70 ELISA pair [C17.8/IL12-IV]

Cat. No. GTX03074**Applications** ELISA, Sandwich ELISA**Reactivity** Mouse**Package**

1 pair

PRODUCT

Product Content

Cat No.	Type	Host	Conjugation	Package
GTX02952-02	Detection antibody	Rat	Biotin	250 µg
GTX02954	Capture antibody	Rat	Unconjugated	250 µg

Applications

Product Note

GTX02954 (C17.8) is able to detect native and recombinant mouse p40 subunit of IL-12 and IL-23. In combination with GTX02952-02 (IL12-IV), it recognizes IL-12 p70 heterodimer (p40 + p35).

Properties

Form Liquid**Storage**


Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

**GTX03074 ELISA Image**

Sandwich ELISA analysis of mouse IL-12 (p70) protein using GTX02954 IL12B / IL12 p40 antibody [C17.8] as coating antibody and GTX02952-02 IL12A / IL12 p35 antibody [IL12-IV] (Biotin) as detecting antibody.



For full product information, images and publications, please visit our [website](#).