

IL12B / IL12 p40 ELISA pair [C15.6/C17.8]

Cat. No. GTX03075

Applications ELISA, Sandwich ELISA

Reactivity Mouse

Package

1 pair

PRODUCT

Product Content

Cat No.	Type	Host	Conjugation	Package
GTX02953	Capture antibody	Rat	Unconjugated	250 µg
GTX02954-02	Detection antibody	Rat	Biotin	250 µg

Applications

Product Note

GTX02953 (C15.6) is able to detect native and recombinant mouse p40 subunit of IL-12 and IL-23. In combination with GTX02954-02 (C17.8), it recognizes p40 monomer, p40 homodimer, IL-12 p70 heterodimer (p40 + p35), and IL-23 heterodimer (p40 + p19).

Properties

Form Liquid

Storage


Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTX03075 ELISA Image

Sandwich ELISA analysis of mouse IL-12/-23 (p40) protein using GTX02953 IL12B / IL12 p40 antibody [C15.6] as coating antibody and GTX02954-02 IL12B / IL12 p40 antibody [C17.8] (Biotin) as detecting antibody.



For full product information, images and publications, please visit our [website](#).