

## IL2 ELISA pair [MT2A91/2C95/MT8G10]

**Cat. No. GTX03081****Applications** ELISA, Sandwich ELISA**Reactivity** Human**Package**

1 pair

## PRODUCT

## Product Content


Cat No.	Type	Host	Conjugation	Package
GTX02975	Capture antibody	Mouse	Unconjugated	250 µg
GTX02976-02	Detection antibody	Mouse	Biotin	250 µg

## Properties

**Form** Liquid**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Note** For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

## DATA IMAGES

**GTX03081 ELISA Image**

Sandwich ELISA analysis of human IL-2 protein using GTX02975 IL2 antibody [MT2A91 + 2C95] as coating antibody and GTX02976-02 IL2 antibody [MT8G10] (Biotin) as detecting antibody.

For full product information, images and publications, please visit our [website](#).