

## IL22 ELISA pair [MT12A3/MT7B27]

**Cat. No. GTX03084****Applications** ELISA, Sandwich ELISA**Reactivity** Human**Package**

1 pair

## PRODUCT

## Product Content


Cat No.	Type	Host	Conjugation	Package
GTX02979	Capture antibody	Rat	Unconjugated	250 µg
GTX02982-02	Detection antibody	Rat	Biotin	250 µg

## Properties

**Form** Liquid**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Note** For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

## DATA IMAGES

**GTX03084 ELISA Image**

Sandwich ELISA analysis of human IL-22 protein using GTX02979 IL22 antibody [MT12A3] as coating antibody and GTX02982-02 IL22 antibody [MT7B27] (Biotin) as detecting antibody.

For full product information, images and publications, please visit our [website](#).