

IL4 ELISA pair [bIL4-I/bIL4-II]

Cat. No. GTX03090

Applications ELISA, Sandwich ELISA

Reactivity Bovine

Package

1 pair

PRODUCT

Product Content

Cat No.	Type	Host	Conjugation	Package
GTX02990	Capture antibody	Mouse	Unconjugated	250 µg
GTX02991-02	Detection antibody	Mouse	Biotin	250 µg

Properties


Form Liquid

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTX03090 ELISA Image

Sandwich ELISA analysis of bovine IL-4 protein using GTX02990 IL4 antibody [bIL4-I] as coating antibody and GTX02991-02 IL4 antibody [bIL4-II] (Biotin) as detecting antibody.

Substrate : TMB



For full product information, images and publications, please visit our [website](#).