

IL4 ELISA pair [11B11/BVD6-24G2]

Cat. No. GTX03092

Applications ELISA, Sandwich ELISA

Reactivity Mouse

Package

1 pair

PRODUCT

Product Content

Cat No.	Type	Host	Conjugation	Package
GTX02989	Capture antibody	Rat	Unconjugated	250 µg
GTX02992-02	Detection antibody	Rat	Biotin	250 µg

Properties


Form Liquid

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTX03092 ELISA Image

Sandwich ELISA analysis of mouse IL-4 protein using GTX02989 IL4 antibody [11B11] as coating antibody and GTX02992-02 IL4 antibody [BVD6-24G2] (Biotin) as detecting antibody.



For full product information, images and publications, please visit our [website](#).