

Interferon alpha ELISA pair [MT1/MT2]

Cat. No. GTX03098

Applications ELISA, Sandwich ELISA

Reactivity Human

Package

1 pair

PRODUCT

Product Content

Cat No.	Type	Host	Conjugation	Package
GTX03014	Capture antibody	Mouse	Unconjugated	250 µg
GTX03016-02	Detection antibody	Mouse	Biotin	250 µg

Applications

Product Note This ELISA pair are able to detect native and recombinant human IFN- α subtypes 2a, 2b, and 2c.

Properties


Form Liquid

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTX03098 ELISA Image

Sandwich ELISA analysis of human IFN- α type 2 protein using GTX03014 Interferon alpha antibody [MT1] as coating antibody and GTX03016-02 Interferon alpha antibody [MT2] (Biotin) as detecting antibody.



For full product information, images and publications, please visit our [website](https://www.genetex.com).