

Interferon gamma ELISA pair [MT17.1/MT307]

Cat. No. GTX03099**Applications** ELISA, Sandwich ELISA**Reactivity** Bovine**Package**

1 pair

PRODUCT

Product Content


Cat No.	Type	Host	Conjugation	Package
GTX02935	Capture antibody	Mouse	Unconjugated	250 µg
GTX02936-02	Detection antibody	Mouse	Biotin	250 µg

Properties

Form Liquid**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Note** For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

**GTX03099 ELISA Image**Sandwich ELISA analysis of bovine IFN- γ protein using GTX02935 Interferon gamma antibody [MT17.1] as coating antibody and GTX02936-02 Interferon gamma antibody [MT307] (Biotin) as detecting antibody.For full product information, images and publications, please visit our [website](#).