

# Interferon gamma ELISA pair [MT13/MT166]

**Cat. No. GTX03100**
**Application** ELISA, Sandwich ELISA

**Reactivity** Dog

**Package**

1 pair

## PRODUCT

### Product Content

Cat No.	Type	Host	Conjugation	Package
GTX02930	Capture antibody	Mouse	Unconjugated	250 µg
GTX02933-02	Detection antibody	Mouse	Biotin	250 µg

## PROPERTIES

**Form** Liquid

**Storage**


Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

## DATA IMAGES



### GTX03100 ELISA Image

Sandwich ELISA analysis of dog IFN-γ protein using GTX02930 Interferon gamma antibody [MT13] as coating antibody and GTX02933-02 Interferon gamma antibody [MT166] (Biotin) as detecting antibody. Substrate : pNPP



For full product information, images and publications, please visit our [website](#).