

Interferon gamma ELISA pair [MT131/MT114]

Cat. No. GTX03101**Application** ELISA, Sandwich ELISA**Reactivity** Cat**Package**

1 pair

PRODUCT

Product Content

Cat No.	Type	Host	Conjugation	Package
GTX02928-02	Detection antibody	Mouse	Biotin	250 µg
GTX02932	Capture antibody	Mouse	Unconjugated	250 µg

PROPERTIES

Form Liquid**Storage**


Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

**GTX03101 ELISA Image**

Sandwich ELISA analysis of cat IFN-γ protein using GTX02932 Interferon gamma antibody [MT131] as coating antibody and GTX02928-02 Interferon gamma antibody [MT114] (Biotin) as detecting antibody. Substrate : pNPP



For full product information, images and publications, please visit our [website](#).