

Interferon gamma ELISA pair [MT166/MT13]

Cat. No. GTX03102

Applications	ELISA, Sandwich ELISA
Reactivity	Horse

Package
1 pair

PRODUCT

Product Content

Cat No.	Type	Host	Conjugation	Package
GTX02930-02	Detection antibody	Mouse	Biotin	250 µg
GTX02933	Capture antibody	Mouse	Unconjugated	250 µg

Properties

Form Liquid

Storage


Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

**GTX03102 ELISA Image**

Sandwich ELISA analysis of horse IFN- γ protein using GTX02933 Interferon gamma antibody [MT166] as coating antibody and GTX02930-02 Interferon gamma antibody [MT13] (Biotin) as detecting antibody.



For full product information, images and publications, please visit our [website](#).