

## Interferon gamma ELISA pair [1-D1K/7-B6-1]

**Cat. No. GTX03104****Applications** ELISA, Sandwich ELISA**Reactivity** Human**Package**

1 pair

## PRODUCT

## Product Content

| Cat No.     | Type               | Host  | Conjugation  | Package |
|-------------|--------------------|-------|--------------|---------|
| GTX02925    | Capture antibody   | Mouse | Unconjugated | 250 µg  |
| GTX02926-02 | Detection antibody | Mouse | Biotin       | 250 µg  |

## Properties

**Form** Liquid**Storage**


Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

## DATA IMAGES

**GTX03104 ELISA Image**Sandwich ELISA analysis of human IFN- $\gamma$  protein using GTX02925 Interferon gamma antibody [1-D1K] as coating antibody and GTX02926-02 Interferon gamma antibody [7-B6-1] (Biotin) as detecting antibody.For full product information, images and publications, please visit our [website](#).