

Interferon gamma ELISA pair [MT126L/7-B6-1]

Cat. No. GTX03105

Applications ELISA, Sandwich ELISA

Reactivity Human, Monkey

Package

1 pair

PRODUCT

Product Content

Cat No.	Type	Host	Conjugation	Package
GTX02926-02	Detection antibody	Mouse	Biotin	250 µg
GTX02929	Capture antibody	Rat	Unconjugated	250 µg

Properties

Form Liquid

Storage


Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTX03105 ELISA Image

Sandwich ELISA analysis of monkey IFN- γ protein using GTX02929 Interferon gamma antibody [MT126L] as coating antibody and GTX02926-02 Interferon gamma antibody [7-B6-1] (Biotin) as detecting antibody.



For full product information, images and publications, please visit our [website](#).