

## TGF beta 1 ELISA pair [MT593/MT517]

## Cat. No. GTX03110

**Applications** ELISA, Sandwich ELISA

**Reactivity** Human

## Package

1 pair

## PRODUCT

## Product Content

| Cat No.     | Type               | Host  | Conjugation  | Package |
|-------------|--------------------|-------|--------------|---------|
| GTX03003-02 | Detection antibody | Mouse | Biotin       | 250 µg  |
| GTX03004    | Capture antibody   | Mouse | Unconjugated | 250 µg  |

## Properties


**Form** Liquid

**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

**Note** For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

## DATA IMAGES



## GTX03110 ELISA Image

Sandwich ELISA analysis of human Latent TGF- $\beta$ 1 protein using GTX03004 TGF beta 1 antibody [MT593] as coating antibody and GTX03003-02 TGF beta 1 antibody [MT517] (Biotin) as detecting antibody.



For full product information, images and publications, please visit our [website](#).