

TNF alpha ELISA pair [MT21A8/MT15B15]

Cat. No. GTX03113

Applications	ELISA, Sandwich ELISA
Reactivity	Human, Monkey

Package
1 pair

PRODUCT

Product Content

Cat No.	Type	Host	Conjugation	Package
GTX03006-02	Detection antibody	Mouse	Biotin	250 µg
GTX03009	Capture antibody	Mouse	Unconjugated	250 µg

Properties

Form Liquid

Storage


Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTX03113 ELISA Image

Sandwich ELISA analysis of monkey TNF-α protein using GTX03009 TNF alpha antibody [MT21A8] as coating antibody and GTX03006-02 TNF alpha antibody [MT15B15] (Biotin) as detecting antibody. Substrate : TMB



For full product information, images and publications, please visit our [website](#).