

## RALBP1 antibody [GT1337]

Cat. No. GTX03249

|              |                   |
|--------------|-------------------|
| Host         | Rabbit            |
| Clonality    | Monoclonal        |
| Isotype      | IgG               |
| Applications | WB, ICC/IF, IHC-P |
| Reactivity   | Human, Mouse, Rat |

Package  
100 µl

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | 1:500 - 1:2000       |
| ICC/IF             | 1:50 - 1:200         |
| IHC-P              | 1:50 - 1:200         |

Not tested in other applications.

|                   |                                                                                 |
|-------------------|---------------------------------------------------------------------------------|
| Calculated MW     | 76 kDa. ( <a href="#">Note</a> )                                                |
| Observed MW (kDa) | 95 kDa. The observed M.W. is based on the following publication. PMID: 16024601 |

## Properties

|               |                                                                                                                                                                                                                            |
|---------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Form          | Liquid                                                                                                                                                                                                                     |
| Buffer        | PBS, 0.05% BSA, 50% Glycerol                                                                                                                                                                                               |
| Preservative  | 0.02% Sodium azide                                                                                                                                                                                                         |
| Storage       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| Concentration | 0.24 mg/ml (Please refer to the vial label for the specific concentration.)                                                                                                                                                |
| Immunogen     | A synthesized peptide derived from human RALBP1.                                                                                                                                                                           |
| Purification  | Purified by affinity chromatography                                                                                                                                                                                        |
| Conjugation   | Unconjugated                                                                                                                                                                                                               |

For full product information, images and publications, please visit our [website](#).


Date 2026 / 10 / 07 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

## Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


## DATA IMAGES



### GTX03249 WB Image

Various whole cell extracts (30 µg) were separated by 7.5% SDS-PAGE, and the membrane was blotted with RALBP1 antibody [GT1337] (GTX03249) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (GTX213110-01) was used to detect the primary antibody.


The observed M.W. is based on the publication: PMID: 16024601



### GTX03249 WB Image

Various whole cell extracts (30 µg) were separated by 7.5% SDS-PAGE, and the membrane was blotted with RALBP1 antibody [GT1337] (GTX03249) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (GTX213110-01) was used to detect the primary antibody.


The observed M.W. is based on the publication: PMID: 16024601



### GTX03249 IHC-P Image

IHC-P analysis of mouse kidney tissue section using GTX03249 RALBP1 antibody [GT1337].

Dilution : 1:100



### GTX03249 WB Image

WB analysis of THP-1 whole cell lysate using GTX03249 RALBP1 antibody [GT1337].

Dilution : 1:1000

Loading : 25µg per lane



For full product information, images and publications, please visit our [website](https://www.genetex.com).