

Mouse IL2 (ECD) protein, mouse IgG2a Fc tag

Cat. No. GTX03268-pro**Species** Mouse**Package**

50 µg

Properties

Form Lyophilized powder**Buffer** Reconstitute with sterile PBS to 100 µg/ml. Lyophilized from 0.2µm-filtered solution in PBS.**Preservative** No preservatives**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Region/Sequence** The extracellular domain of mouse IL-2 (aa 21-169) is fused to the N-terminus of the Fc region of mouse IgG2a.**Expression System** CHO Cells**Purity** > 98% by SDS-PAGE**Endotoxin** < 0.06 EU/µg (determined by LAL assay)**Conjugation** Unconjugated**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).