

Mycobacterium tuberculosis CFP10 + ESAT6 Chimera protein, His tag

Cat. No. GTX03279-pro

Species	Mycobacterium tuberculosis	Package 500 µg
----------------	----------------------------	--------------------------

Properties	
Form	Liquid
Buffer	20mM Phosphate, 0.15mM NaCl
Preservative	No preservatives
Storage	Store as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.
Concentration	2.23 mg/ml (Please refer to the vial label for the specific concentration.)
Expression System	E. coli
Purification	Purified by proprietary method
Purity	95% by SDS-PAGE
Conjugation	Unconjugated
Note	<p>For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.</p> <p>Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.</p>



For full product information, images and publications, please visit our [website](#).