

IDH1 (R132H Mutant) antibody [R132H]

Cat. No. GTX03323

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	IHC-P
Reactivity	Human

Package
500 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	1:25-1:50

Note : Pretreatment : Tris EDTA Buffer pH 9.0 for 30 min at RT

Not tested in other applications.

Product Note

This antibody is specific for R132H mutation of human IDH1 protein and does not cross react with other mutations.

Properties

Form	Liquid
Buffer	Tissue culture supernatant
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Immunogen	KLH-conjugated linear peptide corresponding to human Isocitrate Dehydrogenase 1 (IDH1)
Purification	Unpurified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX03323 IHC-P Image

IHC-P analysis of human anaplastic oligodendroglioma tissue section using GTX03323 IDH1 (R132H Mutant) antibody [R132H].



For full product information, images and publications, please visit our [website](#).