

## Helicobacter pylori antibody

Cat. No. GTX03330

|              |                     |
|--------------|---------------------|
| Host         | Rabbit              |
| Clonality    | Polyclonal          |
| Isotype      | IgG                 |
| Applications | IHC-P               |
| Reactivity   | Helicobacter pylori |

Package  
500 µl

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| IHC-P              | 1:50-1:100           |

**Note : Pretreatment : EDTA Buffer pH 8.0 for 30 min at RT**

Not tested in other applications.

## Product Note

This antibody reacts with the whole H. pylori. The antibody stains the bacteria present on the surface of the epithelium or in the cytoplasm of the epithelial cells.

## Properties

|               |                                                                                         |
|---------------|-----------------------------------------------------------------------------------------|
| Form          | Liquid                                                                                  |
| Buffer        | PBS, 0.1% BSA                                                                           |
| Preservative  | 0.05% Sodium azide                                                                      |
| Storage       | Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. |
| Concentration | Batch dependent (Please refer to the vial label for the specific concentration.)        |
| Immunogen     | Total lysate of H. pylori                                                               |
| Purification  | Purified IgG<br>From serum                                                              |
| Conjugation   | Unconjugated                                                                            |

## Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

## DATA IMAGES



### GTX03330 IHC-P Image

IHC-P analysis of H. pylori infected human tissue section using GTX03330 Helicobacter pylori antibody.



For full product information, images and publications, please visit our [website](#).