

Rickettsia LPS antibody

Cat. No. GTX03396

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Application	ICC/IF, ELISA, IHC
Reactivity	Rickettsia

Package
1 ml

APPLICATION

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	1:10-1:50
ELISA	1:20-1:200
IHC	Assay dependent

Not tested in other applications.

Product Note Reactive with 24 species within the Spotted Fever Group (R. felis weakly). Negative for R. typhi & R. prowazekii.

PROPERTIES

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	100 µg/ml (Please refer to the vial label for the specific concentration.)
Purification	Purified by Protein A or sequential differential precipitations
Conjugation	Unconjugated

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).