

Mouse anti Pig IgA antibody [K60 1F1]

Cat. No. GTX03415

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	ELISA
Reactivity	Pig

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	1:2000-1:200,000

Not tested in other applications.

Product Note The antibody does not cross react with free alpha heavy chain, or with mu or gamma heavy chains.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Porcine IgA
Purification	Protein A purified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX03415 ELISA Image

ELISA analysis of mini pig saliva using serially diluted GTX03415 Mouse anti Pig IgA antibody [K60 1F1].



For full product information, images and publications, please visit our [website](#).