

Human Apolipoprotein CIII protein, His and GST tag

Cat. No. GTX03431-pro

Species Human**Package**

50 µg

Applications


Observed MW (kDa) 34 kDa.

Properties

Form Lyophilized powder**Buffer** Reconstitute with ddH₂O to 0.1-1.0 mg/mL. Do not vortex. Lyophilized from 100mM NaHCO₃ pH8.3, 500mMNaCl, 1mM EDTA, 1mM DTT, 0.01% SKL, 5% Trehalose.**Preservative** ProClin 300**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Region/Sequence** N-terminal His and GST Tagged; Ser21~Ala99**Expression System** E. coli**Purity** >90%**Endotoxin** < 1 EU/µg (Determined by LAL assay)**Conjugation** Unconjugated**Note** For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTX03431-pro Image

SDS-PAGE of GTX03431-pro Human Apolipoprotein CIII protein, His and GST tag.

For full product information, images and publications, please visit our [website](#).