

# Rat lactase protein

**Cat. No. GTX03432-pro**

**Species** Rat

**Package**  
5 test

## PRODUCT

### Summary

It is a purified partial rat intestinal lactase solubilized in SDS-PAGE sample buffer for ready-to-use application. We would recommend using it as GTX47808's positive control.

## Properties

**Form** Liquid

**Buffer** SDS-PAGE sample buffer

**Preservative** No preservatives

**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

**Expression System** Purified from rat small intestinal

**Conjugation** Unconjugated

### Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).