

Human CA125 protein, His Tag

Cat. No. GTX03487-pro

Species Human**Package**

50 µg

Applications


Observed MW (kDa) 26 kDa.

Properties

Form Lyophilized powder**Buffer** Reconstitute with 20mM Tris (pH8.0) and 150mM NaCl to 0.1-1.0 mg/mL. Do not vortex. Lyophilized from 20mM Tris pH8.0, 150mM NaCl, 1mM EDTA, 1mM DTT, 0.01% SKL, 5% Trehalose.**Preservative** Proclin300**Storage** Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.**Region/Sequence** Arg13977~Ser14117 with N-terminal His Tag**Expression System** E. coli**Purity** >95%**Endotoxin** <1.0EU/1µg (determined by the LAL method)**Conjugation** Unconjugated**Note** For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTX03487-pro Image

SDS-PAGE of GTX03487-pro Human CA125 Protein, His Tag.

For full product information, images and publications, please visit our [website](#).