

14-3-3 zeta antibody

Cat. No. GTX03589

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human, Mouse

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:1000
ICC/IF	1:100-1:500
IHC-P	1:50-1:100

Not tested in other applications.

Calculated MW 28 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 150mM NaCl, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Synthesized peptide derived from human 14-3-3 zeta/delta (196-245).
Purification	Purified by antigen-affinity chromatography From serum
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTx03589 ICC/IF Image

ICC/IF analysis of NIH-3T3 cells using GTx03589 14-3-3 zeta antibody. The picture on the right is blocked with the synthesized peptide.



GTx03589 IHC-P Image

IHC-P analysis of human colon carcinoma tissue using GTx03589 14-3-3 zeta antibody. The picture on the right is blocked with the synthesized peptide.



GTx03589 WB Image

WB analysis of K562 cell lysate using GTx03589 14-3-3 zeta antibody. The lane on the right is blocked with the synthesized peptide.



For full product information, images and publications, please visit our [website](#).