

Canine Distemper Virus nucleocapsid antibody [7517]

Cat. No. GTX03593

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG1
Application	ICC/IF, ELISA, Sandwich ELISA
Reactivity	Canine distemper virus

Package
100 µg

APPLICATION

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	1:10-1:50
ELISA	1:20-1:200
Sandwich ELISA	Assay dependent

Note : Capture : GTX03593, Detection : GTX36691

Not tested in other applications.

Product Note These antibodies are non-reactive with canine parvovirus and canine coronavirus.

PROPERTIES

Form	Liquid
Buffer	PBS
Preservative	0.1% sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.1 mg/ml (Please refer to the vial label for the specific concentration.)
Purification	Purified by protein A chromatography or sequential differential precipitations
Purity	>90%
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).



For full product information, images and publications, please visit our [website](#).