

Streptavidin (AP)

Cat. No. GTX03594

Applications IHC-P

Package
100 ml

PRODUCT

Summary

100mL clear cherry red colored solution of streptavidin-conjugated alkaline phosphatase. For immunohistochemical staining of tissue sections. Incubate for 20-30 minutes at room temperature.

Properties

Form Liquid

Buffer 10mM TBS, 0.1% BSA

Preservative 0.1% Proclin

Storage Store at room temperature.

Conjugation Alkaline Phosphatase (AP)

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).