

Tissue Glue

Cat. No. GTX03612

Applications

IHC-P, IHC-Fr

Package

1 kit

PRODUCT

Summary

This product is designed for use with both frozen tissue sections and paraffin-embedded tissue sections for adhering the tissues firmly onto a glass microscope slide.

Properties

Form

Liquid

Storage

Store at 4°C.

Conjugation


Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTX03612 IHC-P Image

IHC-P analysis of human tonsil tissue using GTX01662 Cytokeratin 7 + 8 antibody [CAM5.2] (ready-to-use) with GTX03612 Tissue Glue.



For full product information, images and publications, please visit our [website](#).