

# Tissue Glue

**Cat. No. GTX03612**

## Applications

IHC-P, IHC-Fr

## Package

1 kit

## PRODUCT

### Summary

This product is designed for use with both frozen tissue sections and paraffin-embedded tissue sections for adhering the tissues firmly onto a glass microscope slide.

## Properties

### Form

Liquid

### Storage

Store at 4°C.

### Conjugation


Unconjugated

### Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

## DATA IMAGES



### GTX03612 IHC-P Image

IHC-P analysis of human tonsil tissue using GTX01662 Cytokeratin 7 + 8 antibody [CAM5.2] (ready-to-use) with GTX03612 Tissue Glue.



For full product information, images and publications, please visit our [website](#).