

Purple Hematoxylin

Cat. No. GTX03613
Applications IHC-P

Package
100 ml

PRODUCT

Summary

GTX03613 Purple Hematoxylin is intended for use in the histologic demonstration of nuclear staining. This staining technique is used to make the critical distinction between a normal nucleus and an abnormal one.

Properties

Form Liquid

Storage Store at room temperature.


Conjugation Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTX03613 IHC-P Image

IHC-P analysis of human colon tissue using GTX73737 Epithelial antigen antibody [Ber-EP4] with GTX03613 Purple Hematoxylin.



For full product information, images and publications, please visit our [website](#).

Date 2026 / 01 / 31 Page 1 of 1