

Purple Hematoxylin


Cat. No. GTX03613**Applications** IHC-P**Package**
100 ml**PRODUCT****Summary**

GTX03613 Purple Hematoxylin is intended for use in the histologic demonstration of nuclear staining. This staining technique is used to make the critical distinction between a normal nucleus and an abnormal one.

Properties**Form** Liquid**Storage** Store at room temperature.**Conjugation** Unconjugated**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES**GTX03613 IHC-P Image**

IHC-P analysis of human colon tissue using GTX73737 Epithelial antigen antibody [Ber-EP4] with GTX03613 Purple Hematoxylin.



For full product information, images and publications, please visit our [website](#).