

# Blue Hematoxylin Kit

**Cat. No. GTX03615**

**Applications** IHC-P

**Package**  
100 ml

## PRODUCT

### Summary

Hematoxylin is a histological staining reagent suitable for visualization of nuclei in tissue sections. Blue Hematoxylin, is a blue nuclear counterstain. Bluing Reagent is used following Counterstaining. Both designed for optimal performance when used in immunohistochemistry.

## Properties

**Form** Liquid

**Storage** Store at room temperature.


**Conjugation** Unconjugated

### Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

## DATA IMAGES



### GTX03615 IHC-P Image

IHC-P analysis of human tonsil tissue using Ki67 antibody with GTX03615 Blue Hematoxylin Kit.



For full product information, images and publications, please visit our [website](#).